

Determination of water in vegetable oil using Titrator TitroLine<sup>®</sup> 7500 KF trace

# Application



### Use

The application describes the procedure of the coulometric water determination in vegetable oil products such as olive-, sunflower or similar vegetable oil.

Appliances			
Titrator:	TitroLine 7500 KF trace M1 - 4		
Electrodes			
Generating electrode:	with diaphragm (TZ 1753) or without diaphragm (TZ 1752)		
Reagents			
Use with diaphragm (TZ 1753):			
Anolyt:	Recommended from Sigma Aldrich: 80 ml HYDRANAL-Coulomat AG-H + 20 ml Chloroform. Merck: CombibiCoulomat frit + additional solvent		
Catholyt	Hydranal Coulomat CG for Hydranal reagents; CombiCoulomat frit for Merck		
Additional solvent	The addition of up to 20 % to the anolyt of a long chain alcohol such as decanol or octanol or chloroform is recommended for Merck combicoulomat.		
Use without diaphragm (TZ 1752):			
Reagent:	Recommended is from Sigma Aldrich: Hydranal Coulomat AG-H, from Merck CombiCoulomat fritless		
Additional solvent	The addition of up to 20 % to the reagent of chloroform or a long chain alcohol such as decanol or octanol is recommended.		
Standard	Standards are available from Merck and Sigma Aldrich. Recommended are the standards with concentration of 0.1 %.		



### Description


Set up the unit and fill the reagents as described in the operating manual. Switch on the instrument and wait until the drift is < 10  $\mu$ g/min and stable. For M3 and M4 (generator electrode with diaphragm) it takes sometimes several hours to get a low drift value.

#### Standard and sample Titration

Before you start the sample titration the first time it is recommended to run some tests with a water standard. Standards with certificat in ampoules are recommended instead of pure water.

#### Standard:

- Open the ampoule
- Use a suitable plastic or glass syringe. Depended on the standard use a needle with a diameter between 0.8 mm and 1.2 mm and a length of minimum 70 mm.
- First rinse the syringe 1-2 times with 1 ml each of the standard then draw up slowly the entire ampoule content in the syringe without air-bubbles.
- Place a 100 ml glass beaker (tall form) on a balance, put the syringe inside and weigh it.
- Press tara
- Press the start button on the TL 7500 KF trace
- Inject about 0,75 1,5 ml of the standard into the titration vessel
- Place the syringe inside the glass baker on the balance and read the exact weight from the display/or press the print button for automatic transfer.
- Enter sample ID and sample weight. The titration starts automatically.
- Repeat the determination 2 -3 times.



## Date: 19.12.2013



### Sample:

- Open the sample container
- Use a suitable plastic or glass syringe. Depended on the sample use a needle with a diameter between 0.9 mm and 1.5 mm.
- First rinse the syringe 1-2 times with the sample and then then draw up slowly the sample in the syringe without air-bubbles.
- Place a 100 ml glass beaker (tall form) on a balance, put the syringe inside and weigh it.
- Press tara
- Press the start button on the TL 7500 KF trace
- Inject about 1 2 ml of the sample into the titration vessel
- Place the syringe inside the glass baker on the balance and read the exact weight from the display/or press the print button for automatic transfer.
- Enter sample ID and sample weight. The titration starts automatically.

Note:

If the sample is not homogenous it has to be homogenised before.

# Application



**GLP** documentation Titration graph Olive oil ppm 700,0 600,0 500,0 400,0 300,0 200,0 100,0 0,0 0 s 25,0 50,0 75,0 100,0 125,0 150,0 175,0 200,0 225,0 • ppm/s Method data Method name: Olive oil Titration duration: 3 m 42 s End date: End time: 10:32:59 16.05.13 Titration data Sample ID: Olive oil Weight: 0.68190 g Start drift: 64.7 µg/min End drift: 7.0 µg/min Water: 488.529 µg Water: 716.4 ppm

### Calculation formula

Water:	µg*M*F1/(F2*W)	Mol (M):	1.00000
Factor 1 (F1):	1.0000	Factor 2 (F2):	1.0000
Weight (W):	0.68190 g (m)	Statistics:	Off

# Application



#### Method

#### Method data overall view

Method name: Method type: Olive oil Automatic titration Created at: Last modification: Documentation: 05/16/13 9:59:24 05/16/13 9:59:24 GLP

Start drift: Stop drift (delta):	10.0 µg/min 2.0 µg/min	
Stop drift tolerance:	0.02 ug/min <sup>2</sup>	
Stop delay time:	5 s	
Min. titration time:	60 s	
Max. titration time:	600 s	
Working point:	300 mV	
Control factor:	4	





Hints

If you have any questions concerning the application, you are welcome to contact us.

Literature

优莱博技术 (北京)有限公司

