

Determination of Free and Total SO₂ in Wine and fruit juice according to European regulations

date: 03.04.2013 page 1 from 6



USE

This application note descripe the titration procedure with iodine solution.

APPLIANCES

- Titrator: TL 7000 or TL 7750 M1
- Basic device
- Magnetic stirrer TM 235
- 20 mL Exchange unit WA 20, with amber glass bottle for the titrant, complete
- Option: autosampler TW alpha plus 24 or TW 7400 + 2 x TITRONIC universal

ELECTRODES

Electrode: Pt 1200 with L1NN or Pt 1400 with autosampler

Electrode cable: L 1 NN

REAGENTS

Solvent: water dest.Standardisation: Na2S2O3

Titrant: lodine sol. (I2) 0.025 m or 0.01 m

other reagents
 H2SO4 10 %, NaOH 4 mol/l, KI-solution 5 % and EDTA-Na2

DESCRIPTION

Preparation of Iodine Solution

We recommend ready to use agents. The lodine solutions are also available in ampoules.



A) Determination of the free SO₂

Pipette 50 ml of the sample in a 100 or 150 ml glass beaker (room temperature -> 20 °C), add 3 ml H_2SO_4 10 % + 30 mg EDTA-Na₂, 10 ml KI solution and titrate immediately with the iodine solution.

B) Determination of total SO₂

Add 8 ml 4 m NaOH to the titrated sample A and wait 5 minutes. Add then 10 ml H2SO4 while stirring and titrate immediately with the iodine solution.

Better results could be reached if 20 ml NaOH are added again to the titrated sample and a reaction time of 5 minutes is given. After the addition of H_2SO_4 the titration is started again. Take the summary of both titrations of the total SO2 titration to calculate the result.

Electrode handling

The dead stop electrode can be stored dry after use.

LITERATURE

Amtsblatt der Europäischen Gemeinschaft.


Dr. Alfred Schmitt Aktuelle Weinanalytik Verlag Heller Chemie- und Verwaltungsvorschriften mbH;



Titration Method: SO2 free, page 1

GLP documentation





Method data

Method name: End date:

SO2 in wine 22.02.13

Titration duration:

10:24:18

Titration data

Sample ID: Start µA:

White wine

Pattern:

End time:

25,000 ml

0.192 µA

End µA:

2.793 µA

EP:

1.413 ml/ 2.0 µA

502:

1.8 mg/l

Calculation formula

502:

(EP-B)*T*M*F1/(V*F2)

Mol (M):

32,00000

Blank value (B): Factor 1 (F1):

0.0000 ml 1.0000

Titre (T): Pattern (V):

1.00000000 (m) 25.000 ml (m)

Statistics:

Off

Factor 2 (F2):

1.0000



Titration Method: SO2 free, page 1

Method data overall view

Method name:

Method type: Measured value: Titration mode:

Linear steps:

SO2 in wine Automatic titration

d-stop 0.040 ml Created at:

02/22/13 10:18:41 Last modification: 02/22/13 10:20:53

Documentation:

GLP

Measuring speed / drift:

Initial waiting time: Titration direction: Pretitration:

0 s Increase Off

1 s

Endpoint:

2.0 µA

delta endpoint: Endpoint delay:

1.0 µA 5 s

Polarization voltage:

100 mV

Dosing parameter

Dosing speed: 60.00 % Maximum dosing volume: 20.00 ml

Filling speed:

30 s

Unit values

Unit size: 50ml Unit ID: 10045002 Reagent: Iodid/Iodat Batch ID: no entry Concentration [mol/l]: 1.00000

Determined at: 08/31/12 23:19:38

Expire date; 09/29/12 Opened/compounded: 08/29/12 Test according ISO 8655: 06/01/12

Last modification: 11/22/12 15:47:54



NOTES

If you are using an autosampler then the addition of H_2SO4 and KI have to be done with a piston burette such as TITRONIC universal direct before the titration of the SO_2 .

If you should have any questions concerning application, please contact the Application Department of SI Analytics; tel.: + 49 6131 66 5062 or 5118

优莱博技术(北京)有限公司

地址:北京市朝阳区酒仙桥东路1号院M8号楼C厅3层301室

■E-Mail : info@julabo.cn
Web: www.julabo.cn

