

Potentiometric determination of Silver in Silver or alloys



date: 28.2.2012 Seite 1 von 5



Use

Potentiometric determination of silver in silver or alloys. After dissolving in concentrated nitric acid the sample is titrated with NaCl (or KBr) or HCl using a silver combination electrode.

Appliances

Titrator: TL 6000/7000/7750

- Basic device
- Magnetic stirrer TM 235
- 20 exchange unit WA 20, with brown glass bottle for titrant complete

Electrodes

- Silver combination electrode Ag Cl 62 or AgCl 62 RG
- Electrolyte KNO3 2 mol/l, L 2114 (only for AgCl 62)
- Electrode cable L 1 A

Reagents

- Titrant: NaCl, KCl or HCl 0. Mol/l
- Titer: possible with AgNO3 titrant solution 0.1 mol/l
- Other reagent is nitric acid 65 %

Description

Determination of the exact concentration of the titration (option)

Ready to use titrant solutions are recommended. Also possible is the use of ampules. The titrant solution can be standardised using 10 ml or 15 ml volumes of AgNO3 0.1 mol/l titrant.

For a standard titration pipette 10 - 20 ml of the Ag- standard in a 150 ml beaker, add appr. 80 ml dest water and 1 ml HNO3 (not necessary for HCl titrant) . Place the electrode and burette tip in the sample and start the method. The titration should stops at the equivalence point. The titration consumption (EQ) should be between 10 - 20 ml.

Titration of the sample

Connect the Ag type electrode with the electrode cable at the electrode plug of the TL alpha plus. Do not forget to open the refilling hole of the electrode (only AgCl 62).


Weigh between 180 – 200 mg of the silver sample in a 250 ml glass beaker and add 30 ml of the nitric acid.. Head up the sample with a burner or a heater to dissolve the silver. Dilute the dissolved sample with appr. 100 dest. water. Place the electrode and burette tip in the sample and press the "START" key. The titration stops at the equivalence point.

date: 28.2.2012 Seite 2 von 5



Result example:





Method data

Method name: Silber in % Titration duration: 4 m 34 s
End date: 12.10.12 End time: 16:51:53

Titration data

 Sample ID:
 Charge 3
 Weight:
 0.2018 g

 Start mV:
 348.0 mV
 End mV;
 79.9 mV

EQ: 17.966 ml / 184.6 mV Silber: 96.03 %

Calculation formula

Silber: (EQ1-B)*T*M*F1/(W*F2)

Mol (M): 107.86800

 Blank value (B):
 0.0000 ml
 Titre (T):
 0.10000000

 Factor 1 (F1):
 0.1000
 Weight (W):
 man

 Factor 2 (F2):
 1.0000
 Statistics:
 Off

Device information

Device: TitroLine 7000 Serial number: 10003645

Software version: 1230k Silber_in_%_12_10_12-16_47_19.pdf 1/2

date: 28.2.2012 Seite 3 von 5



Method

The method parameter are very similar to the default method "chloride in %". Only the name and molecular weight have to be changed.

Method data overall view

Method name:Silber in %Created at:10/01/12 16:02:55Method type:Automatic titrationLast modification:10/02/12 18:38:47

Measured value: mV Damping settings: None Titration mode: Dynamic Documentation: GLP

Dynamic: Steep

Measuring speed / drift: User-defined: minimum holding time: 04 s

maximum holding time: 15 s

Measuring time: 03 s

Drift: 10 mV/min

Initial waiting time: 0 s
Titration direction: Decrease
Pretitration: Off
End value: Off
EQ: On (1)

Slope value: Steep Value: 700

Dosing parameter

Dosing speed: 100 % Filling speed; 30 s

Maximum dosing volume: 20.00 ml

Unit values

 Unit size:
 20ml

 Unit ID:
 10039005

 Reagent:
 HCl 0.1 mol/L

 Batch ID:
 no Charge

 Concentration [mol/l]:
 0.10010

Determined at: 12/05/11 19:18:45

Expire date; 08/18/12
Opened/compounded: 09/10/11
Test according ISO 8655: 05/10/11

Last modification: 09/13/12 14:39:31

Device information

Device: TitroLine 7000 Serial number: 10003645

Software version: 1230k Silber_in_%_12_10_12-16_47_19.pdf 2/2

date: 28.2.2012 Seite 4 von 5



Notes

If you have any questions on the application, you can feel free to contact us..

优莱博技术(北京)有限公司

地址:北京市朝阳区酒仙桥东路1号院M8号楼C厅3层301室

■E-Mail : info@julabo.cn
♦Web: www.julabo.cn



date: 28.2.2012 Seite 5 von 5